




## X QC Video (For Semen)

### Product Overview

Push Storage (Capless Design) / USB Interface / Body



### Specification

#### Product Name

LensHooke® X QC Video (For Semen)

#### Size

L57 x W17 x H9 mm

#### Material

Silvery metal (case, plug);  
Black plastic (insulator)

#### Platform

Win7 professional  
Win8.1  
Win10 are recommended.

#### Operation environment

0°C ~ 60°C, 32°F ~ 140°F  
Humidity 20~90%

#### Storage environment



-25°C ~ 85°C, 13°F ~ 185°F  
Humidity 20~90%

#### Transportation environment

-25°C ~ 85°C, 13°F ~ 185°F  
Humidity 20~90%

#### Target Value

Item	Motility (%)	Morphology (%)
V0	0	0
V1	10 ± 2	6 ± 1
V2	40 ± 2	12 ± 1



#### Automatic Rebound Function

it'll automatically rebound while the USB is unplugged.